


Alere INRatio[®]2

PT/INR Monitor



True QC You Can See For Rapid Results That You Can Trust

Only the **Alere INRatio[®]2 PT/INR Monitor** performs normal (QC1) and therapeutic (QC2) assays automatically on every test and provides quantitative, recordable results on the display—true QC you can see, start to finish.

More Than Just A Humidity Check

These two levels of control use a true clotting control measure, following the standards of a split sample laboratory QC system, to provide assurance of testing integrity on every test.

Easy Test Procedure

A simple fingerstick and one unmeasured drop of blood delivers results in about 60 seconds.



Proven Adherence to Clinical Practice Guidelines

The Alere INRatio[®]2 Monitor surpasses the ACCP guidelines for low ISI thromboplastin.³

A Smart Way to Improve Patient Care and Practice Efficiencies

Reduce unnecessary delays

Results in about 60 seconds instead of hours or even days.

Allows prompt medication changes

Counsel your patients face-to-face with dosage changes.

Improve patient satisfaction

Send patients home instead of to a lab for a PT/INR test, or consider patient self testing.

Maximize practice efficiencies

Reduce time-consuming “call-backs”, freeing your staff for more productive duties.

PC Connect[™] Software*

Free Windows software allows you to view, save, print, and export test results.

E-Trainer

Train staff easily with a self-paced, interactive, web or CD based educational tool.

User Friendly Features

New

- Segmented button with audible and tactile feedback
- Print results using the INRatio[®] Printer
- Simple, icon based interface
- Extended battery life supports two hundred tests
- Auto turn-on with strip insertion
- Uses AC or battery power



*PC Connect[™] Software for Professional Use Only



Be empowered to handle the challenges of managing your anticoagulation patients.

The Alere INRatio®2 PT/INR Monitor connects reliable results with practical convenience making it an optimal in-office or home testing solution for anticoagulation patient management.



Alere INRatio®2 Test Strip Advantages



On Board Quality Control

Rigorous quality control checks are performed automatically with every test strip, ensuring a reliable result.



Individually Wrapped Strips, No Refrigeration Required

Individually foil-wrapped strips can be stored at room temperature for up to 12 months.



Small Sample Size

One unmeasured drop of blood applied directly to the strip performs both the test and quality controls.

Alere INRatio®2 PT/INR Professional Testing Kit

CAT 0200431

- 1 Alere INRatio®2 PT/INR Monitor
- 1 User Guide
- 1 User Guide with E-Trainer (CD)
- 2 Quick Reference Guides
- 1 Training Video (DVD)
- 1 Carrying Case
- 10 Disposable Lancets
- 10 Capillary Tubes
- 1 Power Supply Kit

Alere INRatio® PT/INR Test Strips

CAT 0100139

- 48 Test Strips
- 1 Package insert

Reimbursement Information*

CPT Code: 85610QW

Description: Prothrombin Time

(QW modifier indicates a CLIA waived test)

*NOTE: Reimbursement may vary depending on your location/plan. It is the responsibility of the healthcare provider to properly code, document treatment performed, and seek reimbursement.



About Us

By developing new capabilities in near-patient diagnosis, monitoring and health management, Alere gives decision-makers the tools they need to make informed healthcare choices. The result? We're empowering people by putting them on the path to better health.

alere.com | 877-441-7440

References

³ The Seventh ACCP Conference on Anticoagulation Therapy: Evidence-Based Guidelines, Supplement to Chest, Volume 126, No. 3, September 2004: 2085-2095.

©2011 Alere. All rights reserved. The Alere Logo, Alere, INRatio, and PC Connect are trademarks of the Alere group of companies. 1000038vC 10/01

Alere INRatio®2
PT/INR Monitor

